



FLAT SHO ULDER IMPECTOR 19 MM CAP

CAT. NO.	PRODUCT DESCRIPTION	MAT.
1573.062	FLAT SHO ULDER IMPECTOR 19 MM CAP FOR VETERINARY THR	SS

Specification

- Instrument is available in SS.

Indication for use

- The Flat Shoulder Impactor (19 mm Cap) is used in orthopaedic and Total Hip Replacement (THR) surgeries for the controlled impaction of trial or final implants. The 19 mm flat shoulder cap provides stable engagement with the implant surface, allowing accurate seating and alignment of components while minimizing the risk of surface damage or misalignment during insertion.