



IMPACTOR ROD		
CAT. NO.	PRODUCT DESCRIPTION	MAT.
1573.064	CAP IMPECTOR ROD FOR VETERINARY THR	SS
1573.065	ROD IMPECTOR 45 ° FOR VETERINARY THR	SS

Specification

- Instrument is available in SS.

Indication for use

- The Impactor Rod is used in orthopaedic and Total Hip Replacement (THR) surgeries to deliver controlled force for the secure seating of implants or trial components. It provides stability and precision during impaction, ensuring accurate alignment and fixation of prosthetic parts without causing damage to surrounding bone or implant surfaces.

