



IMPECTOR 45 ° POSITION SLEEVE		
CAT. NO.	PRODUCT DESCRIPTION	MAT.
1573.063	IMPECTOR 45* POSITION SLEEVE FOR VETERINARY THR	SS

## Specification

- Instrument is available in SS.

## Indication for use

- The Impactor 45° Position Sleeve is used during Total Hip Replacement (THR) procedures to assist in positioning and impacting the acetabular cup at a precise 45° inclination angle. It ensures accurate alignment and orientation of the implant, promoting optimal joint stability and range of motion while minimizing the risk of improper placement.

