



JIG		
CAT. NO.	PRODUCT DESCRIPTION	MAT.
1573.070	JIG FOR VETERINARY STEM FOR VETERINARY THR	SS

Specification

- Instrument is available in SS.

Indication for use

- The Jig for Veterinary Stem is used during Veterinary Total Hip Replacement (THR) procedures to assist in the accurate alignment and positioning of the femoral stem. It ensures proper orientation and depth during implantation, helping achieve optimal joint stability, biomechanics, and long-term implant performance in veterinary patients.