



NECK		
CAT. NO.	PRODUCT DESCRIPTION	MAT.
800.275S	NECK EXTRA SMALL FOR VETERINARY JOINTS	TT
800.276S	NECK SMALL FOR VETERINARY JOINTS	TT
800.278S	NECK MEDIUM FOR VETERINARY JOINTS	TT
800.280S	NECK LARGE FOR VETERINARY JOINTS	TT
800.282S	NECK EXTRA LARGE FOR VETERINARY JOINTS	TT

## Specification

- Product is available in TT.

## Indication for use

- The neck is used in total hip replacement (THR) or bipolar hip arthroplasty (BHA) to connect the femoral stem to the femoral head, ensuring accurate positioning and smooth joint articulation. It allows for adjustment of offset and leg length, helping to restore normal biomechanics and range of motion of the hip joint.

