



RASP		
CAT. NO.	PRODUCT DESCRIPTION	MAT.
1573.092	RASP EXTRA SMALL FOR VETERINARY THR	SS
1573.094	RASP SMALL FOR VETERINARY THR	SS
1573.096	RASP MEDIUM FOR VETERINARY THR	SS
1573.098	RASP LARGE FOR VETERINARY THR	SS

Specification

- Instrument is available in SS.

Indication for use

- The Rasp is indicated for use in veterinary THR surgeries to accurately contour and enlarge the femoral canal according to the selected implant size. It helps achieve correct alignment, optimal stem fit, and secure seating of the femoral component, ensuring long-term implant stability.