



TRAIL HEAD		
CAT. NO.	PRODUCT DESCRIPTION	MAT.
1573.101	TRAIL HEAD 16 MM EXTRA SMALL	SS
1573.102	TRAIL HEAD 16 MM SMALL	SS
1573.103	TRAIL HEAD 16 MM MEDIUM	SS
1573.104	TRAIL HEAD 16 MM LARGE	SS
1573.105	TRAIL HEAD 16 MM EXTRA LARGE	SS
1573.110	TRAIL HEAD 19 MM SMALL	SS
1573.111	TRAIL HEAD 19 MM MEDIUM	SS
1573.112	TRAIL HEAD 19 MM LARGE	SS
1573.113	TRAIL HEAD 19 MM EXTRA LARGE	SS

Specification

- Instrument is available in SS.

Indication for use

- The Trial Head is used during veterinary Total Hip Replacement procedures to assess the appropriate size and fit of the femoral head component before final implantation. It helps ensure optimal joint stability, range of motion, and articulation with the acetabular cup. The trial head allows surgeons to verify leg length and joint tension intraoperatively, minimizing post-operative complications.

