



TRIAL CUP		
CAT. NO.	PRODUCT DESCRIPTION	MAT.
1573.043	TRAIL CUP 23 MM FOR VETERINARY THR	SS
1573.046	TRAIL CUP 26.5 MM FOR VETERINARY THR	SS
1573.049	TRAIL CUP 29.5 MM FOR VETERINARY THR	SS

Specification

- Instrument is available in SS.

Indication for use

- The Trial Cup is used during veterinary Total Hip Replacement surgery to determine the correct size and positioning of the acetabular component before final implantation. It allows the surgeon to evaluate joint stability, range of motion, and alignment with the femoral head trial. Proper use ensures accurate fit and optimal function of the prosthetic hip joint.