



TRIAL STEM		
CAT. NO.	PRODUCT DESCRIPTION	MAT.
1573.080	TRAIL STEM EXTRA SMALL FOR VETERINARY THR	SS
1573.082	TRAIL STEM SMALL FOR VETERINARY THR	SS
1573.084	TRAIL STEM MEDIUM FOR VETERINARY THR	SS
1573.086	TRAIL STEM LARGE FOR VETERINARY THR	SS

Specification

- Instrument is available in SS.

Indication for use

- The Trial Stem is indicated for temporary intraoperative use in veterinary THR surgeries to evaluate stem size, positioning, and stability before final implantation. It assists in ensuring accurate femoral canal preparation and optimal implant selection, improving the overall fit and function of the prosthesis.