



UNCEMENTED CUP LINNER		
CAT. NO.	PRODUCT DESCRIPTION	MAT.
800.143S	UNCEMENTED CUP LINNER 23.0MM STERILISED FOR VETERINARY	UHMWPE
800.146S	UNCEMENTED CUP LINNER 26.5MM STERILISED FOR VETERINARY	UHMWPE
800.149S	UNCEMENTED CUP LINNER 29.5MM STERILISED FOR VETERINARY	UHMWPE

Specification

- Uncemented Cup Liner is available in UHMWPE material.

Indication for use

- The uncemented cup liner is used in total hip replacement surgeries to provide a low-friction, durable articulation surface between the acetabular cup and femoral head. It ensures smooth joint movement and reduces wear, helping to restore normal hip function. The liner is used in combination with an uncemented acetabular cup in both human and veterinary orthopaedic procedures.