



CAT. NO. SS	Kirschner's Wire (K. Wire)(Sterilised)
241.103S	K. Wire Ø 1.0 mm- L 150mm(Sterilized)
241.121S	K. Wire Ø 1.2 mm- L 150mm(Sterilized)
241.151S	K. Wire Ø 1.5 mm- L 150mm(Sterilized)
241.152S	K. Wire Ø 1.5 mm- L 200mm(Sterilized)
241.181S	K. Wire Ø 1.8 mm- L 150mm(Sterilized)
241.201S	K. Wire Ø 2.0 mm- L 150mm(Sterilized)
241.202S	K. Wire Ø 2.0 mm- L 175mm(Sterilized)
241.203S	K. Wire Ø 2.0 mm- L 200mm(Sterilized)
241.251S	K. Wire Ø 2.5 mm- L 150mm(Sterilized)
241.252S	K. Wire Ø 2.5 mm- L 200mm(Sterilized)
241.301S	K. Wire Ø 3.0 mm- L 150mm(Sterilized)
241.302S	K. Wire Ø 3.0 mm- L 200mm(Sterilized)
241.421S	K. Wire Ø 1.5 mm- L 175mm(Sterilized)
241.426S	K. Wire Ø 2.5 mm- L 175mm(Sterilized)
241.431S	K. Wire Ø 3.0mm - L 175mm(Sterilized)

Specification

- Kirschner's Wire (K.wire)(Sterilised)available in Stainless steel& Titanium.
- Available length150-175 with 1.0,1.2,1.5,1.8,2.0,2.5& 3.0mm dia.

Indication for use

Kirschner's Wire (K.wire)(Sterilised) are used for temporary fixation during some operations.

After definitive fixation they are then removed. The pins are usually removed four to five weeks post operation

