



ANTI TORQUE KEY		
CAT. NO.	PRODUCT DESCRIPTION	MAT.
1573.009	ANTI TORQUE KEY FOR VETERINARY SPINE	SS

Specification

- Product is available in SS.

Indication for use

- The Anti Torque Key is used in spinal and orthopaedic implant procedures to prevent rotational movement of screws or rods during tightening. It provides counter-torque support while the final tightening instrument is engaged. The device ensures stable fixation, accurate implant placement, and reduces the risk of implant displacement during surgery.