



AWL		
CAT. NO.	PRODUCT DESCRIPTION	MAT.
1573.003	AWL FOR VETERINARY SPINE	SS

## Specification

- Product is available in SS.

## Indication for use

- The Awl for Veterinary Spine is used to create an initial pilot hole in the vertebra to facilitate accurate placement of pedicle screws during spinal fixation procedures. It helps ensure correct trajectory and entry point preparation. This instrument supports safe and controlled access to the pedicle. Suitable for spinal stabilization surgeries in small and large animals.