



FINAL TIGHTNER		
CAT. NO.	PRODUCT DESCRIPTION	MAT.
1573.007	FINAL TIGHTNER FOR VETERINARY SPINE	SS

Specification

- Product is available in SS.

Indication for use

- The Final Tightener for Veterinary Spine is used to apply the final locking torque to set screws or locking caps during spinal fixation procedures. It ensures secure and stable fixation of rods within pedicle screws. This instrument provides controlled tightening to prevent loosening and maintain construct stability. Suitable for use in spinal stabilization surgeries for small and large animals.