



### INNI HOLDER

CAT. NO.	PRODUCT DESCRIPTION	MAT.
1573.011	INNI HOLDER FOR VETERINARY SPINE	SS

### Specification

- Product is available in SS.

### Indication for use

- The Inni Holder is used in veterinary spinal fixation procedures to securely hold and guide implants such as screws or rods during insertion. It provides a stable grip that helps maintain accurate alignment in narrow surgical spaces. The instrument enhances control, reduces slippage, and supports precise implant placement throughout the procedure.