



PARALLEL ROD COMPRESSOR		
CAT. NO.	PRODUCT DESCRIPTION	MAT.
1573.008	PARALLEL ROD COMPRESSOR FOR VETERINARY SPINE	SS

Specification

- Product is available in SS.

Indication for use

- The Parallel Rod Compressor is used in veterinary spinal surgery to apply controlled compression between implanted spinal rods. It helps in achieving optimal stabilization and correction of spinal alignment during fixation procedures. The instrument allows precise and parallel force application to support proper fusion and implant positioning. Suitable for use in small and large animal spinal constructs.