



QUICK COUPLING T - HANDLE		
CAT. NO.	PRODUCT DESCRIPTION	MAT.
1573.016	QUICK COUPLING T - HANDLE FOR VETERINARY SPINE	SS

Specification

- Product is available in SS.

Indication for use

- The Quick Coupling T-Handle is used in orthopaedic and spinal procedures to provide controlled manual torque during the insertion or tightening of implants. Its quick-coupling mechanism enables fast and secure attachment of compatible instruments. The T-shaped design offers improved grip and leverage, enhancing precision and reducing hand fatigue during surgical manipulation.

