



ROD HOLDING FORCEPS		
CAT. NO.	PRODUCT DESCRIPTION	MAT.
1573.004	ROD HOLDING FORCEPS FOR VETERINARY SPINE	SS

Specification

- Product is available in SS.

Indication for use

• Rod Holding Forceps are used in spinal fixation procedures to securely grasp, stabilize, and manipulate spinal rods during insertion. They provide a firm, controlled hold to assist in accurate rod placement and alignment. These forceps help maintain rod position while tightening screws or applying corrective maneuvers. Suitable for use in veterinary spinal stabilization surgeries.