



ROD PUSHER

CAT. NO.	PRODUCT DESCRIPTION	MAT.
1573.002	ROD PUSHER FOR VETERINARY SPINE	SS

Specification

- Product is available in SS.

Indication for use

- The Rod Pusher is used during veterinary spinal fixation procedures to accurately guide and position spinal rods into place. It helps maintain alignment while securing rods into pedicle screws. This instrument provides controlled pushing force, ensuring stability during rod insertion. It is suitable for small and large animal spinal corrective or stabilizing surgeries.