



TORQUE QUICK COUPLING HANDLE		
CAT. NO.	PRODUCT DESCRIPTION	MAT.
1573.013	TORQUE QUICK COUPLING HANDLE FOR VETERINARY SPINE	SS

Specification

- Product is available in SS.

Indication for use

- The Torque Quick Coupling Handle is used in orthopaedic and spinal procedures to apply controlled torque during the tightening of surgical implants. Its quick-coupling mechanism allows fast attachment and secure engagement with compatible instruments. The handle ensures consistent torque delivery, improving precision and reducing the risk of overtightening or implant damage.