



| LC-DCP NARROW PLATE 4.5MM | | | |
|---------------------------|-------|-------------|----------------|
| CAT. NO. | HOLES | LENGTH (MM) | THICKNESS (MM) |
| 500.301 | 2 | 34 | 5.2 |
| 500.302 | 3 | 52 | 5.2 |
| 500.303 | 4 | 70 | 5.2 |
| 500.304 | 5 | 88 | 5.2 |
| 500.305 | 6 | 106 | 5.2 |
| 500.306 | 7 | 124 | 5.2 |
| 500.307 | 8 | 142 | 5.2 |
| 500.308 | 9 | 160 | 5.2 |
| 500.309 | 10 | 178 | 5.2 |
| 500.310 | 11 | 196 | 5.2 |
| 500.311 | 12 | 214 | 5.2 |
| 500.312 | 13 | 232 | 5.2 |
| 500.313 | 14 | 250 | 5.2 |
| 500.314 | 15 | 268 | 5.2 |
| 500.315 | 16 | 286 | 5.2 |
| 500.316 | 18 | 323 | 5.2 |
| 500.317 | 20 | 358 | 5.2 |
| 500.318 | 22 | 394 | 5.2 |

Specification

- LC-DCP plate available in both SS and TT materials.
- The plate contains 2-22 holes which has 34-394mm in length and 5.2mm thickness and width 13.5mm. Which is used with 4.5mm cortex screw

Indication for use

- These plates are used for long-bone fractures in small breed dogs and cats.

