



REAMER		
CAT. NO.	PRODUCT DESCRIPTION	MAT.
1573.017	REAMER 17 MM FOR VETERINARY THR	SS
1573.018	REAMER 18 MM FOR VETERINARY THR	SS
1573.019	REAMER 19.5 MM FOR VETERINARY THR	SS
1573.021	REAMER 20.5 MM FOR VETERINARY THR	SS
1573.022	REAMER 22 MM FOR VETERINARY THR	SS
1573.023	REAMER 23.5 MM FOR VETERINARY THR	SS
1573.026	REAMER 26 MM FOR VETERINARY THR	SS
1573.028	REAMER 28 MM FOR VETERINARY THR	SS
1573.030	REAMER 30 MM FOR VETERINARY THR	SS

Specification

- Instrument is available in SS.

Indication for use

- The reamer is used in veterinary Total Hip Replacement (THR) procedures to accurately prepare and shape the acetabular cavity for implant placement. It ensures proper sizing and alignment of the acetabular cup, promoting optimal joint stability and fixation. Suitable for use in small and large animal orthopaedic surgeries.